

Deltalab Sterile Whirl-Pak Sample Bags (2070ml)

Product Code: 200347

Brand: Deltalab

Deltalab's Sterile Whirl-Pak Sample Bags (2070ml) are made of a low density polyethylene blend, and are very resistant and transparent.

Their special polyethylene blend makes unnecessary a higher thickness to achieve higher resistance. They have a leakproof enclosure by several metallic rounded sticks when you turn the closure bands three times.

Used for both liquid and solid samples, these laboratory bags are unique because of their one-piece seam, avoiding the possible risk of the loss of the corners.

They are sterilised by ethylene oxide and are made with materials suitable for alimentary use.

Suitable for liquid nitrogen.

Specifications:

Capacity (ml)	2070
Dimensions (cm)	19 x 38
Thickness (microns)	76
Alimentary Use?	Yes
Case Quantity	500



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com