

Surface Kit

Product Code: 200399

Brand: Deltalab

Surface kits are supplied with a swab and a tube containing a buffered broth with neutralising agents.

Kits are intended for surface sampling for subsequent microbiological analysis.


The kit is a single use product and is sterilised by radiation. The swab is composed of a polystyrene shaft (with breakpoint) and a viscose tip.

The tube incorporates a 'push' cap.

The kit is presented individually wrapped in a peel-pack format.

Specifications:

Shelf Life	18 months
Material	Polystyrene / Viscose
Sterility	Sterile
Quantity	6x 100



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com