

## Deltalab Standard 1.5ml Microtubes (Natural)


**Product Code:** 200400

**Brand:** Deltalab

Made of ultra clear polypropylene, Deltalab's Standard 1.5ml Microtubes (Natural) come with attached flat caps and are autoclavable.

Each microtube can be easily opened and closed with one hand. They are ideal for the storage and freezing of any kind of biological material (cells, sperm, bacteria, etc.) and can be used at temperatures down to  $-100^{\circ}\text{C}$ .

They can also withstand configuration up to 21,000xg equivalent to 15,000rpm.



### Specifications:

<b>Dimensions (Ø x h)</b>	10.86 x 39.06mm
<b>Colour</b>	Natural
<b>Material</b>	Polypropylene
<b>Case Quantity</b>	1000

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)