

Deltalab Graduated Microtubes (1.5ml)

Product Code: 200400P

Brand: Deltalab

Manufactured in polypropylene, Deltalab's Graduated Microtubes (1.5ml) are autoclavable and can be easily opened and closed using their flat caps.

Each microtube also features a frosted area for writing and can withstand centrifugation up to 12,500xg.

Specifications:

Dimensions (Ø x h)	10.2 x 39.9mm
Material	Polypropylene
Case Quantity	12 x 500



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com