

Deltalab L Shaped Spreaders (1 Unit)

Product Code: 200510

Brand: Deltalab

The sterile Deltalab L Shaped Spreaders (1 Unit) are designed for surface spreading.

Made of white polystyrene and available in peel-packs, their base length allows liquid samples to be spread over any type of petri dish and avoids contact with the dish walls.

The curved spreader extremity and rounded corners avoid damaging the media surface while spreading.

The tip bent and rounded edges minimize the possibility of breakage of the solid medium during extension.

It allows their use in places where sterilization before use is not possible (for instance when samples are collected far from the lab).

Specifications:

Presentation	1 unit peel-pack
Total Length (mm)	149
Base Length (mm)	40
Material	Polystyrene
Sterility	Sterile
Quantity	500

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com