

## Deltalab L Shaped Spreaders (I Unit)

**Product Code:** 200510

**Brand:** Deltalab

The sterile Deltalab L Shaped Spreaders (I Unit) are designed for surface spreading.

Made of white polystyrene and available in peel-packs, their base length allows liquid samples to be spread over any type of petri dish and avoids contact with the dish walls.

The curved spreader extremity and rounded corners avoid damaging the media surface while spreading.

The tip bent and rounded edges minimize the possibility of breakage of the solid medium during extension.

It allows their use in places where sterilization before use is not possible (for instance when samples are collected far from the lab).

### Specifications:

<b>Presentation</b>	I unit peel-pack
<b>Total Length (mm)</b>	149
<b>Base Length (mm)</b>	40
<b>Material</b>	Polystyrene
<b>Sterility</b>	Sterile
<b>Quantity</b>	500



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)