

Disposable Sampling Cups

Product Code: 202317

Brand: Deltalab

Disposable polystyrene sample cups, for general sample collection.

Cups are supplied stacked.



Specifications:

Volume (ml)	100
Material	Polystyrene
Mouth Ø (mm)	57.5
Height (mm)	62
Quantity	30x 100

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com