

Yellow 6.7mm Phthalate-Free Anti-Microbial ESD PVC Gas Medical Hose (Vacuum)

Product Code: 208603PF

The Yellow 6.7mm Phthalate-Free Anti-Microbial ESD PVC Gas Medical Hose (Vacuum) provides resistance to bacteria such as MRSA and Klebsiella Pneumoniae, in accordance with BS ISO 22196.

Features

- Contains a low migration plasticiser that meets with the latest chemical toxicity requirements of the European Pharmacopoeia monograph 3.1.2.
- SVHC free
- Will not expose patients to any potentially harmful phthalate vapours
- Does not work harden
- Made from ESD material with a surface resistivity of 1.1×10^{11} ohms/Sq when tested to ASTM D 257
- Manufactured from a wipe down cover containing a slow-release anti-microbial agent that is resistant to many hospital-borne bacteria

Specifications

| | |
|-------------------------------|--------|
| Colour | Yellow |
| Coil Length (m) | 30 |
| Internal Diameter (mm) | 6.7 |
| Gas Type | Vacuum |



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com