

## Blue 6.7mm Phthalate-Free Anti-Microbial ESD PVC Gas Medical Hose (Nitrous Oxide)

**Product Code: 208604PF**

The Blue 6.7mm Phthalate-Free Anti-Microbial ESD PVC Gas Medical Hose (Nitrous Oxide) provides resistance to bacteria such as MRSA and Klebsiella Pneumoniae, in accordance with BS ISO 22196.

### Features

- Contains a low migration plasticiser that meets with the latest chemical toxicity requirements of the European Pharmacopoeia monograph 3.1.2
- SVHC free
- Will not expose patients to any potentially harmful phthalate vapours
- Does not work harden
- Manufactured from a wipe down cover containing a slow-release anti-microbial agent that is resistant to many hospital-borne bacteria

### Specifications

<b>Colour</b>	Blue
<b>Coil Length (m)</b>	30
<b>Gas Type</b>	Nitrous Oxide
<b>Internal Diameter (mm)</b>	6.7

*For more specific information, please consult our sales team.*

*Primehose is manufactured from non-phthalate material; however, it is also manufactured in an environment where phthalates may be present. Therefore, we cannot guarantee that trace elements will not be detectable.*



<b>Material Description</b>
Flexible PVC medical extrusion Gamma Resistant (GIR) compound - non-phthalate plasticised PVC.
<b>Toxicological Information Liner and Cover Material</b>
Free from substances of very high concern as described under REACH regulation. The material complies with the migration requirements of 2002/72/EC and its amendments. Additives used in the compound are listed in EU Regulation 2011/10/EC and in the FDA code of federal regulations Chapter 21.
The cover contains the Biocidal Product Silver Glass Phosphate.
This material is free of Bisphenol A.
<b>Plasticisers &amp; Bonding Agents (Cyclohexanone)</b>
No cyclohexanone is used in the manufacture of these hoses; therefore, there will be no post production leaching of this chemical.
<b>Manufacturing Environment &amp; Storage</b>
Surface contamination will not exceed the requirements stated in BS EN ISO 15001:2001, the product will be fully cling wrapped with the ends left open.

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)