

Medline Sterile Fine Tip Pasteur Pipettes (5ml)

Product Code: 210222

Brand: Medline Scientific

Carry out exact transferring and dispensing of liquids carefully with Medline's Sterile Fine Tip Pasteur Pipettes (5ml).

Made of transparent and inert polyethylene, each pasteur pipette is sterilized by ethylene oxide.



Specifications:

Total Length (mm)	144
Stem Diameter (mm)	5
Drops (ml)	65
Capacity (ml)	5
Material	Polyethylene
Sterility	Sterile
Quantity	500

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com