

## Jeio Tech Funnel Clamp (250ml)

**Product Code:** 23562

**Brand:** Jeio Tech

Mountable, the Jeio Tech Funnel Clamp (250ml) is specially designed for separatory funnels.



### Specifications:

Separatory Funnel Size		250ml
<b>Max. Mountable Funnel Clamps</b>	<b>IST-3075(R)</b>	4
	<b>IST-4075(R)</b>	6
	<b>ISF-7100(R)</b>	15
	<b>ISF-7200(R)</b>	
	<b>ISS-7100(R)</b>	
	<b>ISS-7200(R)</b>	
	<b>ISS-3075(R)</b>	4
	<b>ISS-4075(R)</b>	6

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)