

Jeio Tech Funnel Clamp (500ml)

Product Code: 23563

Brand: Jeio Tech

Mountable, the Jeio Tech Funnel Clamp (500ml) is specially designed for separatory funnels.



Specifications:

Separatory Funnel Size		500ml
Max. Mountable Funnel Clamps	IST-3075(R)	2
	IST-4075(R)	4
	ISF-7100(R)	10
	ISF-7200(R)	
	ISS-7100(R)	
	ISS-7200(R)	
	ISS-3075(R)	3
	ISS-4075(R)	4

Unit 3, Tower Business Park
 Warpsgrove Lane
 Chalgrove
 Oxfordshire
 OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com