

## Jeio Tech Funnel Clamp (2 Litres)

**Product Code:** 23565

**Brand:** Jeio Tech

Mountable, the Jeio Tech Funnel Clamp (2 Litres) is specially designed for separatory funnels.



### Specifications:

Separatory Funnel Size		2 litres
Max. Mountable Funnel Clamps	IST-3075(R)	-
	IST-4075(R)	-
	ISF-7100(R)	5
	ISF-7200(R)	
	ISS-7100(R)	
	ISS-7200(R)	
	ISS-3075(R)	-
	ISS-4075(R)	-

Unit 3, Tower Business Park  
 Warpsgrove Lane  
 Chalgrove  
 Oxfordshire  
 OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)