


## Jeio Tech Test Tube Rack for 4 Tubes (50mm)

**Product Code:** 23593

**Brand:** Jeio Tech

Easy to release and mountable, the Jeio Tech Test Tube Rack for 4 Tubes (50mm) is made of durable stainless steel and provides adjustable angle and tube height support. It can also be used with an optional universal platform.



### Specifications:

<b>No. of Test Tubes</b>		4
<b>Diameter (mm)</b>		50
<b>Mountable Capacity of Racks</b>	<b>IST-3075(R)</b>	2
	<b>IST-4075(R)</b>	4
	<b>ISF-7100(R)</b>	7
	<b>ISF-7200(R)</b>	
	<b>ISS-7100(R)</b>	
	<b>ISS-7200(R)</b>	
	<b>ISS-3075(R)</b>	2
	<b>ISS-4075(R)</b>	4

Unit 3, Tower Business Park  
 Warpsgrove Lane  
 Chalgrove  
 Oxfordshire  
 OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)