

Sterile Pre-Moistened Pointed Foam Swab & Collection Tube

Product Code: 25-88050PF

Brand: Puritan

Pre-Moistened Sterile Collection Sampling Kit for surface sampling

- Intended for environmental surface sampling for the detection of microorganisms
- For sampling surfaces or large equipment parts
- Perfect for spot checks
- Quick-turn-open leak resistant cap
- Swab securely attached to cap
- Available with a pre-moistened round or pointed foam-tipped swab
- Rigid foam paddle tip
- Large tip for maximized collection and elution
- Puritan EnviroMax Plus® brand
- Shelf life - 12 months



Specifications:

Length (mm)	153
Material (Shaft)	Polypropylene
Material (Tip)	Foam
Tip Type (Shape)	Pointed
Sterility	Sterile
Quantity	250

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com