

Deltalab Sterile Water Sampling Bottles with Sodium Thiosulfate (500ml) - Packs of 24

Product Code: 282320

Brand: Deltalab

Deltalab's Sterile Water Sampling Bottles with Sodium Thiosulfate (500ml) - Packs of 24 are specifically designed for examining water that is used for drinking and swimming pools, along with other types of water where chlorine may alter the arrangement of the sample while being transported.

Each water sampling bottle is labelled specifying description, code, lot, and expiry date.

Features:

- Minor sides are gripped for a better handling.
- Larger sides are flat for labelling.
- Tamper evident red cap with an internal joint.
- Hermetical closure.

Specifications:

Capacity (ml)	500
Dimensions (mm)	83 x 65 x 135
Mouth Diameter (mm)	28
Unitary Weight (empty bottle)	40g
Material	Polyethylene
Sterility	Sterile
Quantity	24



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com