

Deltalab Sterile Water Sampling Bottles with Sodium Thiosulfate (1000ml)

Product Code: 282324.BU

Brand: Deltalab

Deltalab's 1000ml Water Sampling Bottles are specifically designed for examining water that is used for drinking, swimming pools, along with other types of water whereby chlorine may alter the arrangement of the sample whilst being transported.

Each water sampling bottle is sterile and has a printed label specifying description, code, lot number, and expiry date.

Individually wrapped

Features:

- Minor sides are gripped for a better handling.
- Larger sides are flat for labelling.
- Tamper evident red cap with an internal joint.
- Hermetical closure.

Specifications:

Capacity (ml)	1000
Dimensions (mm)	83 x 65 x 235
Mouth Diameter (mm)	28
Unitary Weight (empty bottle)	61g
Material	Polyethylene
Sterility	Sterile
Quantity	68



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com