

## Deltalab Sterile Water Sampling Bottles with Sodium Thiosulfate (1000ml)

**Product Code:** 282324

**Brand:** Deltalab

Deltalab's 1000ml Water Sampling Bottles are specifically designed for examining water that is used for drinking, swimming pools, along with other types of water whereby chlorine may alter the arrangement of the sample whilst being transported.

Each water sampling bottle is sterile and has a printed label specifying description, code, lot number, and expiry date.

### Features:

- Minor sides are gripped for a better handling.
- Larger sides are flat for labelling.
- Tamper evident red cap with an internal joint.
- Hermetical closure.

### Specifications:

<b>Capacity (ml)</b>	1000
<b>Dimensions (mm)</b>	83 x 65 x 235
<b>Mouth Diameter (mm)</b>	28
<b>Unitary Weight (empty bottle)</b>	61g
<b>Material</b>	Polyethylene
<b>Sterility</b>	Sterile
<b>Quantity</b>	68



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)