

Deltalab Sterile Water Sampling Bottles Without Sodium Thiosulfate (500ml) - Packs of 24

Product Code: 282330

Brand: Deltalab

Deltalab's Sterile Water Sampling Bottles Without Sodium Thiosulfate (500ml) - Packs of 24 feature tamper-evident red screw caps, which assures the bottle's sterility before use. These laboratory containers can be used to sample water for physico-chemical analysis, along with other sampling requiring sterile bottles.

Features:

- Minor sides are gripped for a better handling.
- Larger sides are flat for labelling.
- Tamper evident red cap with an internal joint.
- Hermetical closure.

Specifications:

Capacity (ml)	500
Dimensions (mm)	83 x 65 x 135
Mouth Diameter (mm)	28
Unitary Weight (empty bottle)	40g
Material	Polyethylene
Sterility	Sterile
Quantity	24

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com