

## Deltalab Sterile Water Sampling Bottles Without Sodium Thiosulfate (500ml) - Packs of 111

**Product Code:** 282333

**Brand:** Deltalab

Deltalab's Sterile Water Sampling Bottles Without Sodium Thiosulfate (500ml) - Packs of 111 feature tamper-evident red screw caps, which assures the bottle's sterility before use. These laboratory containers can be used to sample water for physico-chemical analysis, along with other sampling requiring sterile bottles.

### Features:

- Minor sides are gripped for a better handling.
- Larger sides are flat for labelling.
- Tamper evident red cap with an internal joint.
- Hermetical closure.

### Specifications:

<b>Size (ml)</b>	500
<b>Dimensions (mm)</b>	83 x 65 x 135
<b>Mouth Diameter (mm)</b>	28
<b>Unitary Weight (empty bottle)</b>	40g
<b>Material</b>	Polyethylene
<b>Sterility</b>	Sterilised by radiation
<b>Quantity</b>	111



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)