

Deltalab Sterile Water Sampling Bottles Without Sodium Thiosulfate (1000ml) - Packs of 68

Product Code: 282334

Brand: Deltalab

Deltalab's Sterile Water Sampling Bottles Without Sodium Thiosulfate (1000ml) - Packs of 68 feature tamper-evident red screw caps, which assures the bottle's sterility before use. These laboratory containers can be used to sample water for physico-chemical analysis, along with other sampling requiring sterile bottles.

Features:

- Minor sides are gripped for a better handling.
- Larger sides are flat for labelling.
- Tamper evident red cap with an internal joint.
- Hermetical closure.

Specifications:

Size (ml)	1000
Dimensions (mm)	83 x 65 x 235
Mouth Diameter (mm)	28
Unitary Weight (empty bottle)	61g
Material	Polyethylene
Sterility	Sterilised by radiation
Quantity	68



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com