

## Deltalab Sterile Water Sampling Jars - Without Sodium Thiosulfate (1000ml)

**Product Code:** 282351

**Brand:** Deltalab

Deltalab's Sterile Water Sampling Jars - Without Sodium Thiosulfate (1000ml) can be used to sample water for physico-chemical analysis as well as for other sampling requiring sterile bottles.

### Features:

- Transparent. Made of PET (Polyethylene Terephthalate).
- Wide mouth and grips at two sides, easing the sample collection.
- Red cap made of polyethylene with an internal liner made of plexan with a label that performs as tamper evident.

### Specifications:

<b>Internal Mouth Diameter (mm)</b>	55
<b>Capacity (ml)</b>	1000
<b>Jar Weight (g)</b>	65
<b>Sterility</b>	Sterilised by radiation
<b>Material(s)</b>	Polyethylene Terephthalate (jar) and Polyethylene (cap)
<b>Quantity</b>	48

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)