

## Deltalab Tamper Evident Polyethylene Bottles (125ml)

**Product Code:** 292303

**Brand:** Deltalab

The translucent cylindrical Deltalab Tamper Evident Polyethylene Bottles (125ml) are made of high density polyethylene, and come with a PP screw cap and LDPE insert plug. They are specially designed for the storage of alcohol.

Supplied uncapped, the tamper evident caps are ribbed for a better handling while working with gloves. Insert plugs are dispensed pierced. The bottles feature no graduations.



### Specifications:

Capacity (ml)	125
Neck Diameter (mm)	19
Body Diameter (mm)	45
Height (mm)	106
Unitary Weight (g)	10
Material	Polyethylene
Case Quantity	140

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)