

Deltalab Plain Snappable Polystyrene & Viscose Sterile Swabs (Peel-Packed)


Product Code: 300202

Brand: Deltalab

Sterilized by ethylene oxide, Deltalab's Plain Snappable Polystyrene & Viscose Sterile Swabs (Peel-Packed) can be used when samples do not need to be transported.

Specifications:

Material(s)	Snappable Polystyrene & Viscose
Packaging	Peel-Packed
Sterility	Sterile
Quantity	1000



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com