

Deltalab Plain Non-Sterile Swabs (Wooden Shaft Only)

Product Code: 300205


Brand: Deltalab

Suitable for gynaecology when taking endocervical cells while using the speculum, Deltalab's Non-Sterile Swabs (Wooden Shaft Only) have straight-cut splinter free ends.

Supplied in bulk packs.

Specifications

Shaft Dimensions (mm)	147 ± 0.3 x 2.2
Tip Diameter (mm)	± 0.3
Sterility	Non-Sterile
Quantity	1000



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com