

## Deltalab Wooden & Cotton Sterile Swabs (2 Units)

**Product Code:** 300210

**Brand:** Deltalab

Sterilized by ethylene oxide, Deltalab's Wooden & Cotton Sterile Swabs (2 Units) consist of two sterile swabs wrapped in peel-packs.

One swab is designed for cleaning the sampling area. The other swab is designed for sample collection.

For use when samples do not need to be transported.

### Specifications:

<b>Packaging</b>	Peel-Pack
<b>Sterility</b>	Sterile
<b>Quantity</b>	1000

