

Deltalab Plain Non-Sterile Swabs (Aluminium Shaft & Cotton Tip)


Product Code: 300243

Brand: Deltalab

Deltalab's Non-Sterile Swabs (Aluminium Shaft & Cotton Tip) are suitable for gynaecology when taking endocervical cells while using the speculum.

Specifications

Shaft Dimensions (mm)	150 x 0.9
Tip Dimensions (mm)	± 1.5
Sterility	Non-Sterile
Quantity	1000



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com