

Deltalab Plain Sterile Swabs in Round Tubes (Aluminium Shaft and Viscose Tip)

Product Code: 300253

Brand: Deltalab

Sterilised by ethylene oxide, Deltalab's Plain Sterile Swabs in Round Tubes (Aluminium Shaft and Viscose Tip) are dry swabs equipped with shockproof round bottom polypropylene tubes. Each tube has a label to securely seal the cap.

Specifications:

Tube Dimensions (mm)	13 x 165
Material(s)	Aluminium (shaft) Viscose (tip) Polypropylene (tube)
Sterility	Sterile
Quantity	500

