

Deltalab Sterile Cary Blair Swabs (Wood & Cotton)

Product Code: 300280

Brand: Deltalab

Deltalab's Sterile Cary Blair Swabs (Wood & Cotton) are a modification of Stuart media and are specially formulated for faecal collection and anaerobic transport.

All swabs are sterilised by radiation.
Peel-packed.

Specifications

Peel-Pack Dimensions (mm)	38 x 210
Material(s)	Wood and Cotton
Expiry Date	30 months from sterilisation date
Case Weight	8.50
Case Volume	0.056
Sterility	Sterile
Quantity	600

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com