

Deltalab Sterile Amies Liquid Swabs (Twisted Aluminium & Viscose)


Product Code: 300281T

Brand: Deltalab

Deltalab's Sterile Amies Liquid Swabs (Twisted Aluminium & Viscose) are a variation of the traditional Amies media swabs where the media is presented on liquid form (in a sponge) without agar.

These laboratory swabs are suitable for the conversation of pathogens such as *Trichomonas Vaginalis*, *Streptococcus Pyogenes*, and *Shigella Flexneri*.

All swabs are sterilised by radiation.
Peel-packed.



Specifications

Peel-Pack Dimensions (mm)	38 x 210
Material(s)	Twisted Aluminium and Viscose
Expiry Date	30 months after sterilisation date
Quantity	540

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com