


## Deltalab Sterile Amies Swabs (Snappable Polystyrene & Viscose)

**Product Code:** 300287

**Brand:** Deltalab

Deltalab's Sterile Amies Swabs (Snappable Polystyrene & Viscose) are a modification of Cary Blair media, which is at the same time a modification of Stuart media. It assures the viability of organisms such as *Trichomonas Vaginalis*, *Streptococcus Pyogenes*, and *Streptococcus Pneumoniae*.

Swabs are sterilised by radiation.  
Peel-packed.



### Specifications

<b>Peel-Pack Dimensions (mm)</b>	38 x 210
<b>Material(s)</b>	Snappable Polystyrene and Viscose
<b>Expiry Date</b>	30 months after sterilisation date
<b>Quantity</b>	100

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)