

## Deltalab Sterile Stuart Swabs (Aluminium & Cotton)

**Product Code:** 300291

**Brand:** Deltalab

Deltalab's Sterile Stuart Swabs (Aluminium & Cotton) are modified Stuart media transport swabs that allow the conservation and transportation of a large number of pathological microorganisms (e.g. neisseria gonorrhoeae, haemophilus influenzae).

The most unstable organisms will remain visible for up to 24 hours and others for several days. The media is reduced due to the presence of thioglycolate, which difficults the enzymatic reactions of the bacteria.

The multiplication of the bacteria is prevented due to the lack of nitrogen in the media.

### Specifications:

<b>Peel-Pack Dimensions (mm)</b>	38 x 210
<b>Material(s)</b>	Aluminium and Cotton
<b>Expiry Date</b>	30 months from sterilisation date
<b>Sterility</b>	Sterile
<b>Quantity</b>	600



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)