

Deltalab Sterile Stuart Swabs (Snappable Polystyrene & Viscose)

Product Code: 300295

Brand: Deltalab

Deltalab's Sterile Stuart Swabs (Snappable Polystyrene & Viscose) are modified Stuart media transport swabs that allow the conservation and transportation of a large number of pathological microorganisms (e.g. neisseria gonorrhoeae, haemophilus influenzae).

The most unstable organisms will remain visible for up to 24 hours and others for several days. The media is reduced due to the presence of thioglycolate, which diffcults the enzymatic reactions of the bacteria.

The multiplication of the bacteria is prevented due to the lack of nitrogen in the media.

Specifications:

Peel-Pack Dimensions (mm)	38 x 210
Material(s)	Snappable Polystyrene and Viscose
Expiry Date	30 months from sterilisation date
Sterility	Sterile
Quantity	600

