

Deltalab Sterile Viral Swabs (Snappable Polystyrene & Polyester)

Product Code: 300297

Brand: Deltalab

Deltalab's Sterile Viral Swabs preserve specimens (e.g. Papiloma, Pseudorrabia) during transportation to the laboratory.

The media maintains the sample viable for up to 72 hours and in some cases for a longer period of time.

The addition of antimicrobial substances inhibit the growth of bacteria and fungi. Ideal for nasal, pharyngeal, ocular, skin and mucus samples.

Media is supplied in its liquid form, in a sponge, suitable for direct extensions on slides.

All swabs are sterilised by radiation.

Specifications:

Material(s)	Snappable Polystyrene and Viscose
Expiry Date	18 months from sterilisation date
Case Weight	7.75
Case Volume	0.057
Sterility	Sterile
Quantity	100

