

## Deltalab Sterile Chlamydia Swabs (Snappable Polystyrene & Polyester)


**Product Code:** 300299

**Brand:** Deltalab

Deltalab's Sterile Chlamydia Swabs (Snappable Polystyrene & Polyester) are liquid media transport swabs suitable for detecting chlamydia.

For cervical samples, it is recommended to firstly use a dry swab to clear the cervical channel.

The media is supplied in liquid form, in a sponge, suitable for direct extensions on slides. Sterilised by radiation.



### Specifications:

<b>Material(s)</b>	Snappable Polystyrene and Polyester
<b>Expiry Date</b>	12 months from sterilisation
<b>Case Weight</b>	7.80
<b>Case Volume</b>	0.056
<b>Sterility</b>	Sterile
<b>Quantity</b>	100

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)