

Deltalab 4ml Round Bottom Polystyrene Tubes (11 x 70mm)

Product Code: 300300

Brand: Deltalab


Manufactured from high technology, Deltalab's 4ml Round Bottom Polystyrene Tubes (11 x 70mm) feature new generation moulds to ensure reproducibility.

Each test tube is made in one piece to ensure uniformity and dimensional accuracy from tube to tube.

Suitable for centrifugation up to 7500xg, all tubes can be used in various laboratories.

Specifications:

Dimensions (mm)	11 x 70
Volume (ml)	4
Internal Diameter (mm)	9.1
External Diameter (mm)	10.8
Total Length (mm)	69.4
Maximum Volume (ml)	4
Material	Polystyrene
Quantity	3000



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com