

Deltalab 5ml Round Bottom Polystyrene Tubes (12 x 75mm - Thick Wall)

Product Code: 300800.2

Brand: Deltalab

Manufactured from high technology, Deltalab's 5ml Round Bottom Polystyrene Tubes (12 x 75mm - Thick Wall) feature new generation moulds to ensure reproducibility.

Each test tube is made in one piece to ensure uniformity and dimensional accuracy from tube to tube.

Suitable for centrifugation up to 7500xg, all tubes can be used in various laboratories.

Specifications:

Dimensions (mm)	12 x 75 - thick wall
Volume (ml)	5
Internal Diameter (mm)	10
External Diameter (mm)	11.8
Total Length (mm)	74.5
Maximum Volume (ml)	5.3
Material	Polystyrene
Quantity	4000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com