

## Deltalab 5ml Graduated Round Bottom Polystyrene Tubes (12 x 75mm)

**Product Code:** 300800

**Brand:** Deltalab

Manufactured from high technology, Deltalab's 5ml Graduated Round Bottom Polystyrene Tubes (12 x 75mm) feature new generation moulds to ensure reproducibility.


Each test tube is made in one piece to ensure uniformity and dimensional accuracy from tube to tube.

Suitable for centrifugation up to 7500xg, all tubes can be used in various laboratories.

These tubes are graduated at 1, 2.5, and 4ml.

### Specifications:

<b>Dimensions (mm)</b>	12 x 75
<b>Volume (ml)</b>	5
<b>Internal Diameter (mm)</b>	10.3
<b>External Diameter (mm)</b>	11.8
<b>Total Length (mm)</b>	74.5
<b>Maximum Volume (ml)</b>	5.5
<b>Material</b>	Polystyrene
<b>Quantity</b>	4000



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)