

## Deltalab Sterile Polystyrene Culture Tubes (17 x 100mm)

**Product Code:** 300808

**Brand:** Deltalab

Deltalab's Sterile Polystyrene Culture Tubes (17 x 100mm) are supplied with either a two position ribbed polyethylene cap that can be left loose work for aerobic work or seal for anaerobic cultures.

They are biologically inert, exempt from mould release agents, and withstand up to 1400xg and 70°C.

All tubes are packaged in self-standing resealable zip-lock bags.



### Specifications:

<b>Dimensions (mm)</b>	17 x 100
<b>Material</b>	Polystyrene
<b>Case Quantity</b>	1000

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)