

Deltalab 11ml Capped and Labelled Polystyrene Tubes (16 x 100mm)

Product Code: 300908

Brand: Deltalab

Suitable for centrifugation up to 3000xg, Deltalab's 11ml Capped and Labelled Polystyrene Tubes (16 x 100mm) can be used for various laboratory applications (e.g. cell culture, cell separation).

Specifications:

| | |
|-------------------------------------|----------|
| Tube Dimensions (mm) | 16 x 100 |
| Height - Tube & Cap (mm) | 110.4 |
| Volume (ml) | 11 |
| Quantity | 2000 |



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com