

Deltalab 11ml Capped and Labelled Polystyrene Tubes (16 x 100mm)

Product Code: 300908

Brand: Deltalab

Suitable for centrifugation up to 3000xg, Deltalab's 11ml Capped and Labelled Polystyrene Tubes (16 x 100mm) can be used for various laboratory applications (e.g. cell culture, cell separation).

Specifications:

Tube Dimensions (mm)	16 x 100
Height - Tube & Cap (mm)	110.4
Volume (ml)	11
Quantity	2000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com