

Deltalab Sterile Eppendorf Pipette Tip Rack (0.5 - 10µl)

Product Code: 301-09

Brand: Deltalab

Deltalab's Sterile Pipette Tip Rack (0.5 - 10µl) has a hinged lid and blue base. It can hold up to 96 tips.

It is made of autoclavable polypropylene up to 121°C.

Specifications:

Description	96 tips rack
Material	Polypropylene
Sterility	Sterile
Quantity	10



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com