

Deltalab 12ml Conical Polystyrene Tubes (16 x 102mm) - No Graduations

Product Code: 301202

Brand: Deltalab

Deltalab's 12ml Conical Polystyrene Tubes (16 x 100mm) can be used for cell centrifugation and other cell culture applications.

Each conical tube has an external rim with no internal graduations.

Specifications:

Autoclavable	No
Dimensions (mm)	16 x 102
Material	Polystyrene
Volume (ml)	12
Internal Diameter (mm)	14.5
External Diameter (mm)	19.1/16.4
Total Length (mm)	100.8
Max. Volume (ml)	12.9
Quantity	2500

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com