

Deltalab 12ml Conical Polystyrene Tubes (17 x 105mm) - Max. Volume 16.4ml

Product Code: 301212

Brand: Deltalab

Deltalab's 12ml Conical Polystyrene Tubes (17 x 105mm) - Max. Volume 16.4ml can be used for cell centrifugation and other cell culture applications.

Each conical tube is flared at the top for easy pouring and has a special sedimentation base.

Specifications:

Autoclavable	No
Dimensions (mm)	17 x 105
Material	Polystyrene
Volume (ml)	12
Internal Diameter (mm)	15.4/19.2
External Diameter (mm)	17.1/22
Total Length (mm)	105
Max. Volume (ml)	16.4
Quantity	1500



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com