


Deltalab EUROTUBO 12ml Round Bottom Screw Cap Tubes (Polystyrene)

Product Code: 301402

Brand: Deltalab

Made of transparent polystyrene, Deltalab's EUROTUBO 12ml Round Bottom Screw Tubes (Polystyrene) feature a green cap made of high density polyethylene. The external skirt allows the tubes to remain free-standing. Each tube is supplied screwed.



Specifications:

Description	Polystyrene
Recommended Volume (ml)	12
Dimensions (mm)	15 x 102
Autoclavable	No
Sterility	Non-Sterile
Max. Volume (ml)	14.4
External Cap (Ø mm)	20.9
External Tube (Ø mm)	16.3
Internal Tube (Ø mm)	14.4
Length with Cap (Ø mm)	103.9
Length without Cap (Ø mm)	102.5
Quantity	1500

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com