

Deltalab 2mm Non-Calibrated Nichrome Wire Loops

Product Code: 302762

Brand: Deltalab


Deltalab's 2mm Non-Calibrated Nichrome Wire Loops are used for microbiology.

Comparing them to the traditional loops made of platinum and rhodium, these loops have a totally smooth surface and provide rapid cooling after heat sterilization.

Each loop comes without a holder, saving money and space. They are supplied in groups of 25 units in one tube, specifying the code, description, and lot.

Specifications:

Description	Non-calibrated loop 2mm Ø
Dimensions (mm)	75
Quantity	25



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com