

Medline Standard Short and Disposable Loops (1 µl)

Product Code: 302773

Brand: Medline Scientific


Manufactured in flexible polystyrene, Medline's Standard Short and Disposable Loops (1 µl) are designed smooth and flexible in order to facilitate effortless streaking without damaging the gel surface.

Each loop has a rounded tip and are all colour coded. Sterilised by ethylene oxide, they are disposable which eliminates the risk of cross contamination. They are also free of any lubricants and do not contaminate any media.

They are supplied in convenient sterile peel pouch bags with the batch and expiry date clearly displayed.

Specifications:

Style	1 µl loop (rigid type); sterile
Height (mm)	200
Weight (g)	29
Colour	Blue
Material	Polystyrene
Packaging	20 ud./peel-pack, 50 peels/box
Quantity	1000



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com