

Deltalab 4.5ml Special UV Grade Spectrophotometer Cuvettes (4 Sides Transparent)

Product Code: 303100

Brand: Deltalab


Made of UV grade PMMA for accuracy throughout UV/Vis range (280 to 800), Deltalab's 4.5ml Special UV Grade Spectrophotometer Cuvettes (4 Sides Transparent) are disposable and suitable for most open spectrophotometers. They carry out homogeneous measures, specially of the surface crossed by the light beam, assuring an optimum transmission level on the whole visible spectral.

The material used avoids any possible measurement interference. Due to the strict quality control during the manufacture process, a high reliability is assured. The maximum absorption variations are $\pm 1\%$.

The two sides not crossed by the light beam are ribbed to an easy identification of the cuvette position inside the spectrophotometer measurement chamber, resulting in easy positioning and removal.

They are supplied in a dust proof, expandable polystyrene box (100 units per box) with a lid.

With four clear sides, these cuvettes are ideal for fluorometry and nephelometry.



Specifications:

Description	4.5ml 4 sides transparent
Dimensions (mm)	12.55 x 12.65 x 44.55 (± 0.1 mm)
Light Path (mm)	10
Quantity	500

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com