

Deltalab Mauve Polyethylene Caps for 15 - 16mm Tubes (with Internal Chamber)

Product Code: 305008


Brand: Deltalab

Ideal for capping and recapping 15 - 16mm blood collection tubes, Deltalab's Mauve Polyethylene Caps for 15 - 16mm Tubes (with Internal Chamber) are designed to be placed and removed securely. Each cap has an internal chamber for more liquid filling, thus compensating for the loss of space inside the tube caused by the insertion of a cap. They are ergonomical and have a ribbed surface. Their shape is slightly conical for better handling.

These leakproof caps also feature an external flap to avoid contact between liquid and fingers when opening the tube, prevent dripping when the cap is left on the bench top, and provide easy single-handed opening.

Specifications:

Colour	Mauve
Internal Diameter (mm)	13.95
External Diameter (mm)	19.30
Total Length (mm)	13.15
Material	Polyethylene
Sterility	Non-Sterile
Quantity	2500



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com