


Deltalab Natural Polyethylene Caps for 11 - 12mm Tubes (with Internal Chamber)

Product Code: 305302

Brand: Deltalab

Ideal for capping and recapping 11 - 12mm blood collection tubes, Deltalab's Natural Polyethylene Caps for 11 - 12mm Tubes (with Internal Chamber) are designed to be placed and removed securely. Each cap has a small internal chamber and two internal ribs to clamp firmly on the tube. They are also light-weighted and ergonomical.



Specifications:

Colour	Natural
Internal Diameter (mm)	10
External Diameter (mm)	11.85
Total Length (mm)	11.30
Material	Polyethylene
Sterility	Non-Sterile
Quantity	2000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com